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Emphases for proper use:

1. **It is required to replace the product filters (square compressor filter and round bacterial filter) every 3000 hours. Filters replacement is a condition for warranty realization. In the absence of such replacement, warranty would not be granted.**
2. The system is an oxygen concentrator of **1-5 liters per minute**, which means that **93%±3%** oxygen will be achievable at all flow levels between **1 to 5 liters per minute**. Passing the **5-liter per minute** followed by non-linear drop in oxygen level. Chaban's recommendation is not to operate the system above **4.5 liters per minute**.
3. **The system flow is represented by the upper part of the ball in the flow meter**. When the upper part of the ball touches the 5-liter line, the concentrator output is **5 liters**.
4. Chaban's recommendation is to use distilled water in the humidifier bottle. Cleaning the humidifier with vinegar for scale removal shall be perform **once a month**. A blocked humidifier bottle is one of the main causes of system malfunctions and low flow rate. The scale blocks the humidifier bottle inner filter and therefore interfere the oxygen flow rate.
*****Note: Chaban recommend humidifier bottle replacement every 6 months.**
5. The green connector supplied with the product serves as a backup connector in the absence of humidifier bottle or as an alternative if the patient prefer dry oxygen.
The green connector should be kept!

*****Note: The above document is not to replace the user manual!**

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