

V-type Animal Operating Table

TTVOT-01



Product Description

1. The height of the operating table is lifted and lowered by electric control, and the lower tank is equipped in the middle.
2. This working table can be tilted by -5° ~ 15° to the left and right, and mechanically operated manually.
3. The front and back of the platform can be tilted by 45° respectively, with electronic operation.
4. The whole machine is compact, reliable and reasonable, and easy to operate.
5. The countertop is made of 304 stainless steel, which is resistant to high temperature, corrosion and rust.
6. The working surface temperature adjustment range is between 0-60 degrees, and the automatic thermostatic function can be adjusted as required temperature.
7. Built-in heating protector to ensure the 100% safety.
8. Designed four fastener around the table, easy for fixing and reliable.



Product parameters (units are mm)

Dimensions(L x W x H): 1400*650*760-1060mm



Material Description

The main body of the operating table is made of 304 stainless steel material, with stainless steel detachable mesh, convenient and practical High-intensity lifter, The national standard stainless steel 304 material is anti-corrosion, acid-proof and rust-free, with electric lifter, the performance is more stable.

